

Product datasheet for **HP213135**

CBWD1 Human qPCR Primer Pair (NM_018491)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55871
Forward Sequence:	GAAACGGTTGCCTCTGCTGTTC
Reverse Sequence:	AGGGTCTGCTAATCCAGTGGTC
ACCN:	NM_018491 , NM_018491.1 , NM_018491.2 , NM_018491.3 , BC009573 , BC005996 , BC013432 , BC037212 , BC132759 , BU621970 , NM_018491.5
UniProt ID:	Q9BRT8
Synonyms:	COBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)