

Product datasheet for **HP213161**

CSGALNACT2 Human qPCR Primer Pair (NM_018590)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55454
Forward Sequence:	TCATCTCACAGTGGTGTATTTTGG
Reverse Sequence:	GCACCCACATTTAGTCCTCGTC
ACCN:	NM_018590 , NM_018590.1 , NM_018590.2 , NM_018590.3 , NM_018590.4 , BC030268 , BX102005 , BX640967 , BX641073 , BX647369 , NM_018590.5
UniProt ID:	Q8N6G5
Synonyms:	beta4GalNAcT; ChGn-2; CHGN2; GALNACT-2; GALNACT2; PRO0082
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)