

Product datasheet for **HP213172**

TROY (TNFRSF19) Human qPCR Primer Pair (NM_018647)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55504
Forward Sequence:	GGTGCATTCTGCAGCCAGTCTT
Reverse Sequence:	CAGGCATCTGAAAACCTGCCAC
ACCN:	<u>NM_018647</u> , <u>NM_018647.1</u> , <u>NM_018647.2</u> , <u>NM_018647.3</u> , <u>BC047321</u> , <u>BC047321.1</u> , <u>BC035796</u> , <u>BF968836</u> , <u>BF971156</u> , <u>BG674480</u> , <u>BI494770</u> , <u>BM801444</u> , <u>NM_018647.5</u>
UniProt ID:	<u>Q9NS68</u>
Synonyms:	TAJ; TAJ-alpha; TRADE; TROY
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)