

Product datasheet for **HP213173**

NOP10 Human qPCR Primer Pair (NM_018648)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55505
Forward Sequence:	GATGGGACAACAGACCTGCTCA
Reverse Sequence:	TTAAGGGACCCTCAGAGGACAG
ACCN:	BC063023 , NM_018648 , NM_018648.1 , NM_018648.2 , NM_018648.3 , BC063023.1 , BC008886 , BM554889 , BQ218536 , NM_018648.4
UniProt ID:	Q9NPE3
Synonyms:	DKCB1; NOLA3; NOP10P
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)