

Product datasheet for **HP213191**

NSMase2 (SMPD3) Human qPCR Primer Pair (NM_018667)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55512
Forward Sequence:	GAAGCACACCTCAGGACCAAAG
Reverse Sequence:	CAGCCAGTCCTGAAGCAGGTC
ACCN:	<u>NM_018667</u> , <u>NM_018667.1</u> , <u>NM_018667.2</u> , <u>NM_018667.3</u> , <u>BC012727</u> , <u>BC041828</u> , <u>BC112238</u> , <u>BC143631</u> , <u>BC143634</u> , <u>BI196095</u>
UniProt ID:	<u>Q9NY59</u>
Synonyms:	NSMASE2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)