

Product datasheet for **HP213196**

ACCN4 (ASIC4) Human qPCR Primer Pair (NM_018674)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55515
Forward Sequence:	GGAGATCATGCTGGACATCCAG
Reverse Sequence:	GCTGTGGATCTGCACCCGAATA
ACCN:	NM_018674 , NM_018674.1 , NM_018674.2 , NM_018674.3 , NM_018674.4 , NM_018674.5 , BC010439 , BC010439.1 , BC031812 , BC073912
UniProt ID:	Q96FT7
Synonyms:	ACCN4; BNAC4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)