

## Product datasheet for **HP213211**

### PRDM5 Human qPCR Primer Pair (NM\_018699)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 11107  |
| Forward Sequence: | ACTCTGAGGAGAGACCGTTCCA   |
| Reverse Sequence: | AGCATCGCAGTGATGGCACTTG   |
| ACCN:             | <u><a href="#">NM_018699</a></u> , <u><a href="#">NM_018699.1</a></u> , <u><a href="#">NM_018699.2</a></u> , <u><a href="#">NM_018699.3</a></u> , <u><a href="#">BC047115</a></u> , <u><a href="#">BC066942</a></u> , <u><a href="#">BC121037</a></u> , <u><a href="#">BC121038</a></u> , <u><a href="#">BM563688</a></u> , <u><a href="#">BX646150</a></u> , <u><a href="#">NM_018699.4</a></u> |
| UniProt ID:       | <u><a href="#">Q9NQX1</a></u>  |
| Synonyms:         | BCS2; PFM2   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.   |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)