

## Product datasheet for **HP213221**

### RCC2 Human qPCR Primer Pair (NM\_018715)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55920
Forward Sequence:	CCAAACGTGGTTGTACGAGACG
Reverse Sequence:	GGACCATCTCATCCTTCTGCTC
ACCN:	<a href="#">NM_018715</a> , <a href="#">NM_018715.1</a> , <a href="#">NM_018715.2</a> , <a href="#">NM_018715.3</a> , <a href="#">BC053908</a> , <a href="#">BC053908.1</a> , <a href="#">BC004933</a> , <a href="#">BC042141</a> , <a href="#">BC142946</a> , <a href="#">BQ774471</a> , <a href="#">NM_018715.4</a>
UniProt ID:	<a href="#">Q9P258</a>
Synonyms:	TD-60
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)