

Product datasheet for **HP213313**

SLC38A2 Human qPCR Primer Pair (NM_018976)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54407
Forward Sequence:	AGCCAACAGCTCTGTACCTGC
Reverse Sequence:	GGAAGAACAGCAGGATGACAGAC
ACCN:	NM_018976 , NM_018976.1 , NM_018976.2 , NM_018976.3 , NM_018976.4 , BC040342 , BC040342.1 , BC016909 , BC029379 , NM_018976.5
UniProt ID:	Q96QD8
Synonyms:	ATA2; PRO1068; SAT2; SNAT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)