

Product datasheet for **HP213316**

WNK1 Human qPCR Primer Pair (NM_018979)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	65125
Forward Sequence:	CAGTGTGGAACCAGAGGGTGAT
Reverse Sequence:	TTTGCTGTGGCATGGAAGACGC
ACCN:	BC021121 , NM_018979 , NM_018979.1 , NM_018979.2 , NM_018979.3 , BC130467 , BC013629 , BC035146 , BC044600 , BC071959 , BC094862 , BC130469 , BC141881 , BC172444 , BJ997509 , BM716053 , BX644268 , NM_018979.4
UniProt ID:	Q9H4A3
Synonyms:	HSAN2; HSN2; KDP; p65; PPP1R167; PRKWNK1; PSK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)