

Product datasheet for **HP213347**

CSGLCAT (CHPF2) Human qPCR Primer Pair (NM_019015)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54480
Forward Sequence:	AGTTGGAGCGGGCTTACAGTGA
Reverse Sequence:	CAGCACCTCAAAGCGAGAGTGT
ACCN:	<u>NM_019015</u> , <u>NM_019015.1</u> , <u>NM_019015.2</u> , <u>BC052944</u> , <u>BC064509</u>
UniProt ID:	<u>Q9P2E5</u>
Synonyms:	chPF-2; ChSy-3; CSGLCA-T; CSGLcAT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)