

Product datasheet for **HP213358**

DHX29 Human qPCR Primer Pair (NM_019030)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54505
Forward Sequence:	GATGTCCTGGTAGTATGCCCTAC
Reverse Sequence:	CCAGCCAAACATCTCGGTAAGTG
ACCN:	NM_019030 , NM_019030.1 , NM_019030.2 , NM_019030.3 , BC056219 , BC056219.1 , BC035988 , BX648101 , BX648269 , BX648688 , BX649135 , NM_019030.4
UniProt ID:	Q7Z478
Synonyms:	DDX29
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)