

Product datasheet for **HP213362**

PCDH18 Human qPCR Primer Pair (NM_019035)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54510
Forward Sequence:	GGGCAATATCAGCCAAGACCAAG
Reverse Sequence:	ACTGTCTCCTGCCTCACTGTCA
ACCN:	<u>NM_019035</u> , <u>NM_019035.1</u> , <u>NM_019035.2</u> , <u>NM_019035.3</u> , <u>NM_019035.4</u> , <u>BC093815</u> , <u>BC022227</u> , <u>BC071581</u> , <u>BC143361</u> , <u>BG541979</u> , <u>BQ017212</u> , <u>BX097655</u> , <u>BX647786</u> , <u>NM_019035.5</u>
UniProt ID:	<u>Q9HCL0</u>
Synonyms:	PCDH68L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)