

## Product datasheet for **HP213370**

### WDR44 Human qPCR Primer Pair (NM\_019045)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54521
Forward Sequence:	TTGGCAAGTGCTACCAGTGGAG
Reverse Sequence:	CTCCCAGAGTTGGATCTAGGTC
ACCN:	<a href="#">NM_019045</a> , <a href="#">NM_019045.1</a> , <a href="#">NM_019045.3</a> , <a href="#">NM_019045.4</a> , <a href="#">BC028697</a> , <a href="#">BC028697.2</a> , <a href="#">BC013924</a> , <a href="#">BX470418</a> , <a href="#">BX537578</a> , <a href="#">BX538087</a> , <a href="#">BX648150</a> , <a href="#">NM_019045.5</a>
UniProt ID:	<a href="#">Q5JSH3</a>
Synonyms:	RAB11BP; RPH11; SYM-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)