

Product datasheet for **HP213378**

FAM193B Human qPCR Primer Pair (NM_019057)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54540
Forward Sequence:	TCAACAGCACGGAGCCCAAAGT
Reverse Sequence:	GGCTTCCAGAAACACGATTTGCC
ACCN:	NM_019057 , NM_019057.3 , NM_019057.2 , BC008784 , BC071696 , BC086859 , BF768924 , BQ228026 , BX431446 , BX440663 , BX647587
UniProt ID:	Q96PV7
Synonyms:	DKFZp3131142; FLJ10404; KIAA1931
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)