

Product datasheet for **HP213383**

EML4 Human qPCR Primer Pair (NM_019063)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27436
Forward Sequence:	TGCAACTGGACAGATAGCTGGC
Reverse Sequence:	AGGCATCCTACTCCACGCTCAA
ACCN:	<u>NM_019063</u> , <u>NM_019063.1</u> , <u>NM_019063.2</u> , <u>NM_019063.3</u> , <u>NM_019063.4</u> , <u>BC146799</u> , <u>BC008685</u> , <u>BC023522</u> , <u>BC038351</u> , <u>BC067777</u> , <u>BC104647</u> , <u>BC140845</u> , <u>BX647693</u> , <u>BX647731</u>
UniProt ID:	<u>Q9HC35</u>
Synonyms:	C2orf2; ELP120; EMAP-4; EMAPL4; ROPP120
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)