

Product datasheet for **HP213394**

UGT1A10 Human qPCR Primer Pair (NM_019075)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54575
Forward Sequence:	GGGCAAAATCCCTCAGACAGTC
Reverse Sequence:	AGCATGGGTGATAAAGGCACGG
ACCN:	<u>NM_019075</u> , <u>NM_019075.1</u> , <u>NM_019075.2</u> , <u>BC020971</u> , <u>BC020971.1</u> , <u>BC053576</u> , <u>BC069210</u>
UniProt ID:	<u>Q9HAW8</u>
Synonyms:	GNT1; hUG-BR1; UDPGT; UGT-1A; UGT-1J; UGT1; UGT1-01; UGT1-10; UGT1.1; UGT1.10; UGT1A; UGT1A1; UGT1J
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)