

## Product datasheet for **HP213434**

### HOXD8 Human qPCR Primer Pair (NM\_019558)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3234
Forward Sequence:	GTCCAGTGGTAATATTGGCGAGG
Reverse Sequence:	GCGACTGTAGTTTGTCTTCCTC
ACCN:	<a href="#">NM_019558</a> , <a href="#">NM_019558.1</a> , <a href="#">NM_019558.2</a> , <a href="#">NM_019558.3</a> , <a href="#">BC038709</a> , <a href="#">BC038709.1</a> , <a href="#">BC090853</a> , <a href="#">BI764317</a> , <a href="#">BQ429720</a> , <a href="#">NM_019558.4</a>
UniProt ID:	<a href="#">P13378</a>
Synonyms:	HOX4; HOX4E; HOX5.4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)