

Product datasheet for **HP213437**

RNF20 Human qPCR Primer Pair (NM_019592)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56254
Forward Sequence:	GCTGGAGAACTTCGGCAAGAC
Reverse Sequence:	CAAGACGGAGAACTGTGACTGC
ACCN:	<u>NM_019592</u> , <u>NM_019592.1</u> , <u>NM_019592.2</u> , <u>NM_019592.4</u> , <u>NM_019592.5</u> , <u>NM_019592.6</u> , <u>BC110584</u> , <u>BC016691</u> , <u>BC032618</u> , <u>BC063115</u> , <u>BC110585</u> , <u>BC152309</u> , <u>NM_019592.7</u>
UniProt ID:	<u>Q5VTR2</u>
Synonyms:	BRE1; BRE1A; hBRE1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)