

Product datasheet for **HP213438**

KIAA1434 (GPCPD1) Human qPCR Primer Pair (NM_019593)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56261
Forward Sequence:	GGAGTGTGGTTATGGCTTGACAG
Reverse Sequence:	CCTGAACACTACGTGCTCACTGAG
ACCN:	<u>NM_019593</u> , <u>NM_019593.1</u> , <u>NM_019593.2</u> , <u>NM_019593.3</u> , <u>BC027588</u> , <u>BC027588.2</u> , <u>BG674365</u> , <u>BI561344</u> , <u>BP421674</u> , <u>NM_019593.5</u>
UniProt ID:	<u>Q9NPB8</u>
Synonyms:	EDI3; GDE5; GDPD6; PREI4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)