

Product datasheet for **HP213439**

LRR8A Human qPCR Primer Pair (NM_019594)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56262
Forward Sequence:	ACCTCCAGAACCTAGCCATCAC
Reverse Sequence:	CAGTGACTGCAGCACGTTGTTG
ACCN:	<u>NM_019594</u> , <u>NM_019594.1</u> , <u>NM_019594.2</u> , <u>NM_019594.3</u> , <u>BC051322</u> , <u>BC051322.1</u> , <u>BC000775</u> , <u>BX376144</u>
UniProt ID:	<u>Q8IWT6</u>
Synonyms:	AGM5; HsLRR8A; LRR8; SWELL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)