

## Product datasheet for **HP213460**

### EIF4ENIF1 Human qPCR Primer Pair (NM\_019843)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56478
Forward Sequence:	GTGAGCACAATGAAGGCAAGTGG
Reverse Sequence:	GCAGGATGTTCTTAGGGACAGG
ACCN:	<u><a href="#">NM_019843</a></u> , <u><a href="#">NM_019843.1</a></u> , <u><a href="#">NM_019843.2</a></u> , <u><a href="#">NM_019843.3</a></u> , <u><a href="#">BC033028</a></u> , <u><a href="#">BC033028.1</a></u> , <u><a href="#">BC032941</a></u> , <u><a href="#">BC012734</a></u> , <u><a href="#">NM_019843.4</a></u>
UniProt ID:	<u><a href="#">Q9NRA8</a></u>
Synonyms:	4E-T; Clast4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)