

## Product datasheet for **HP213463**

### CCL28 Human qPCR Primer Pair (NM\_019846)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56477
Forward Sequence:	AGCTGTTGCACGGAGGTTTCAC
Reverse Sequence:	ACAGCAGCCAAGTCACAATCCC
ACCN:	<a href="#">NM_019846</a> , <a href="#">BC062668</a> , <a href="#">BC062668.1</a> , <a href="#">BC069532</a> , <a href="#">BX486632</a>
UniProt ID:	<a href="#">Q9NRJ3</a>
Synonyms:	CCK1; chemokine (C-C motif) ligand 28; MEC; MEC, CCK1, SCYA28, MGC71902; Scya28; small inducible cytokine A28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)