

Product datasheet for **HP213475**

DIABLO Human qPCR Primer Pair (NM_019887)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 56616 |
| Forward Sequence: | GGAGGAAGATGAAGTGTGGCAG |
| Reverse Sequence: | GCCATCTCTGAAAGACCAACTGC |
| ACCN: | NM_019887 , NM_019887.1 , NM_019887.2 , NM_019887.3 , NM_019887.4 , NM_019887.5 , BC011909 , BC011909.2 , BC004417 , BC046209 , BC095468 , BF313712 , BF724951 , BM746689 , NM_019887.6 |
| UniProt ID: | Q9NR28 |
| Synonyms: | DFNA64; SMAC |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)