

Product datasheet for **HP213478**

INPP5E Human qPCR Primer Pair (NM_019892)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56623
Forward Sequence:	GACCTCATCTGGTCTGCTCAG
Reverse Sequence:	CTGAGGTGAAGTGGGACGTGAT
ACCN:	<u>NM_019892</u> , <u>NM_019892.1</u> , <u>NM_019892.2</u> , <u>NM_019892.3</u> , <u>NM_019892.4</u> , <u>NM_019892.5</u> , <u>BC028032</u> , <u>BC028032.1</u> , <u>BC110356</u> , <u>BE044535</u> , <u>BE314630</u> , <u>NM_019892.6</u>
UniProt ID:	<u>Q9NRR6</u>
Synonyms:	CORS1; CPD4; JBTS1; MORMS; pharbin; PPI5PIV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)