

Product datasheet for **HP213485**

SLC2A4RG Human qPCR Primer Pair (NM_020062)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56731
Forward Sequence:	TCCAGTGTCTGTGGAAGAGCTG
Reverse Sequence:	CCACACCAACATCCAGCTCTGT
ACCN:	NM_020062 , NM_020062.1 , NM_020062.2 , NM_020062.3 , BC052306 , BC052306.1 , BC001402 , BC017446 , BC028349 , BM983611 , NM_020062.4
UniProt ID:	Q9NR83
Synonyms:	GEF; HDBP-1; HDBP1; Si-1-2; Si-1-2-19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)