

Product datasheet for **HP213496**

SPHK2 Human qPCR Primer Pair (NM_020126)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56848
Forward Sequence:	GAGGAAGCTGTGAAGATGCCTG
Reverse Sequence:	GAGCAGTTGAGCAACAGGTCGA
ACCN:	NM_020126 , NM_020126.1 , NM_020126.2 , NM_020126.3 , NM_020126.4 , BC006161 , BC006161.1 , BC010671 , BI224962 , NM_020126.5
UniProt ID:	Q9NRA0
Synonyms:	SK-2; SK 2; SPK-2; SPK 2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)