

Product datasheet for **HP213499**

TCIM Human qPCR Primer Pair (NM_020130)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56892
Forward Sequence:	AGCCTCTCGTAAGAAAGCCGTG
Reverse Sequence:	CTTGGCTCTCTCCTGCTTTC
ACCN:	NM_020130 , NM_020130.1 , NM_020130.2 , NM_020130.3 , NM_020130.4 , BC021672 , BC021672.1 , BC020623 , BM976768 , BX434108
UniProt ID:	Q9NR00
Synonyms:	C8orf4; TC-1; TC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)