

Product datasheet for **HP213504**

BDH2 Human qPCR Primer Pair (NM_020139)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56898
Forward Sequence:	CACAACCAAGGCAGCCGTGATT
Reverse Sequence:	CTTGTAGAGATGGCGTATCAACTG
ACCN:	<u>NM_020139</u> , <u>NM_020139.1</u> , <u>NM_020139.2</u> , <u>NM_020139.3</u> , <u>BC095414</u> , <u>BC095414.1</u> , <u>BC001953</u> , <u>BC037277</u> , <u>BI544872</u> , <u>NM_020139.4</u>
UniProt ID:	<u>Q9BUT1</u>
Synonyms:	DHRS6; EFA6R; PRO20933; SDR15C1; UCPA-OR; UNQ6308
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)