

Product datasheet for **HP213528**

NMUR2 Human qPCR Primer Pair (NM_020167)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56923
Forward Sequence:	CTGCTGTGTTCAACCTCGTCCA
Reverse Sequence:	ATGCTGGGAGTGCCACTGTTTG
ACCN:	BC016938 , NM_020167 , NM_020167.1 , NM_020167.2 , NM_020167.3 , NM_020167.4 , BC016938.2 , BC067776
UniProt ID:	Q9GZQ4
Synonyms:	FM-4; FM4; NMU-R2; NMU2R; TGR-1; TGR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)