

## Product datasheet for **HP213548**

### MFF Human qPCR Primer Pair (NM\_020194)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56947
Forward Sequence:	CAAGGTTCCAGGCACCGATTTC
Reverse Sequence:	GCGACAAAATGCCACGAGCAGA
ACCN:	<u><a href="#">NM_020194</a></u> , <u><a href="#">NM_020194.1</a></u> , <u><a href="#">NM_020194.2</a></u> , <u><a href="#">NM_020194.3</a></u> , <u><a href="#">NM_020194.4</a></u> , <u><a href="#">NM_020194.5</a></u> , <u><a href="#">BC000797</a></u> , <u><a href="#">BC093024</a></u> , <u><a href="#">BI550249</a></u> , <u><a href="#">BQ653061</a></u>
UniProt ID:	<u><a href="#">Q9GZY8</a></u>
Synonyms:	C2orf33; EMPF2; GL004
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)