

Product datasheet for **HP213550**

XAB2 Human qPCR Primer Pair (NM_020196)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56949
Forward Sequence:	CGACACGAGAACTACGATGAGG
Reverse Sequence:	AGCATGGACCAGACCTTCAGTG
ACCN:	<u>NM_020196, NM_020196.1, NM_020196.2, BC007208, BC007208.1, BC008778, NM_020196.3</u>
UniProt ID:	<u>Q9HCS7</u>
Synonyms:	HCNP; HCRN; NTC90; SYF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)