

Product datasheet for **HP213554**

PRTFDC1 Human qPCR Primer Pair (NM_020200)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56952
Forward Sequence:	GGACCATGAAAGCACTACTCAGC
Reverse Sequence:	GGTCTAAAGCCGTCACCTTCTGG
ACCN:	<u>NM_020200</u> , <u>NM_020200.1</u> , <u>NM_020200.2</u> , <u>NM_020200.3</u> , <u>NM_020200.4</u> , <u>NM_020200.5</u> , <u>NM_020200.6</u> , <u>BC008662</u> , <u>BC008662.2</u> , <u>BQ003226</u> , <u>BX471710</u> , <u>NM_020200.7</u>
UniProt ID:	<u>Q9NRG1</u>
Synonyms:	HHGP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)