

## Product datasheet for **HP213561**

### SEMA4B Human qPCR Primer Pair (NM\_020210)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10509
Forward Sequence:	ACCGTGAGACACAGCAGTGGTA
Reverse Sequence:	GTGGTCCTTGAGGAAGTTCAGC
ACCN:	<a href="#">NM_020210</a> , <a href="#">NM_020210.1</a> , <a href="#">NM_020210.2</a> , <a href="#">NM_020210.3</a> , <a href="#">BC110853</a> , <a href="#">BC110853.1</a> , <a href="#">BC010701</a> , <a href="#">BC017658</a> , <a href="#">BC019013</a> , <a href="#">BC094849</a> , <a href="#">BC126438</a> , <a href="#">BC144141</a> , <a href="#">BX537671</a> , <a href="#">NM_020210.4</a>
UniProt ID:	<a href="#">Q9NPR2</a>
Synonyms:	SEMAC; SemC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)