

## Product datasheet for **HP213575**

### PPAN Human qPCR Primer Pair (NM\_020230)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56342
Forward Sequence:	CGACCACAACATCACAGAGCTG
Reverse Sequence:	TCCTGGACCTTGATGAGCTGCA
ACCN:	<a href="#">NM_020230</a> , <a href="#">NM_020230.1</a> , <a href="#">NM_020230.2</a> , <a href="#">NM_020230.3</a> , <a href="#">NM_020230.4</a> , <a href="#">NM_020230.5</a> , <a href="#">NM_020230.6</a> , <a href="#">BC033202</a> , <a href="#">BC033202.1</a> , <a href="#">BC000535</a> , <a href="#">BC009833</a> , <a href="#">NM_020230.7</a>
UniProt ID:	<a href="#">Q9NQ55</a>
Synonyms:	BXDC3; SSF; SSF-1; SSF1; SSF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)