

Product datasheet for **HP213576**

KTEL1 (POGLUT1) Human qPCR Primer Pair (NM_020231)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56983
Forward Sequence:	CGCCAGAAGGAGTCAGGTTCAA
Reverse Sequence:	CTGGAGATGCCTCCTCGGAAAG
ACCN:	NM_020231 , NM_020231.3 , BC030614 , BC030614.2 , BC048810 , BP347179
UniProt ID:	Q8NBL1
Synonyms:	CLP46, MDSRP, C3orf9, MDS010, hCLP46, KDELCL1, MGC32995
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)