

## Product datasheet for **HP213584**

### KIF15 Human qPCR Primer Pair (NM\_020242)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56992
Forward Sequence:	TTGGACCAAACAGCAGGAAGAGC
Reverse Sequence:	GACTACTCGCAGGTCATGTACC
ACCN:	<u><a href="#">NM_020242</a></u> , <u><a href="#">NM_020242.1</a></u> , <u><a href="#">NM_020242.2</a></u> , <u><a href="#">BC117174</a></u> , <u><a href="#">BC143752</a></u> , <u><a href="#">NM_020242.3</a></u>
UniProt ID:	<u><a href="#">Q9NS87</a></u>
Synonyms:	HKLP2; KLP2; KNSL7; NY-BR-62
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)