

Product datasheet for **HP213594**

VGLuT1 (SLC17A7) Human qPCR Primer Pair (NM_020309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57030
Forward Sequence:	GCAAGTACATCGAGGACGCCAT
Reverse Sequence:	GCCACGATGATGGCATAGACTG
ACCN:	NM_020309 , NM_020309.1 , NM_020309.2 , NM_020309.3 , BC059379 , BC059379.1 , BG701511 , BQ636968 , BU732143 , NM_020309.4
UniProt ID:	Q9P2U7
Synonyms:	BNPI; VGLUT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)