

Product datasheet for **HP213595**

MAX binding protein (MNT) Human qPCR Primer Pair (NM_020310)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4335
Forward Sequence:	ACGTACTGGAGATTGACCGCGT
Reverse Sequence:	CCTCGTCTATGTTGTCCTCACC
ACCN:	<u>NM_020310</u> , <u>NM_020310.1</u> , <u>NM_020310.2</u> , <u>BC117563</u> , <u>BC117563.1</u> , <u>BC033186</u> , <u>BM469815</u> , <u>BM820333</u> , <u>BM876168</u> , <u>BM876428</u> , <u>BQ675978</u> , <u>NM_020310.3</u>
UniProt ID:	<u>Q99583</u>
Synonyms:	bHLHd3; lncRNA-HAL; MAD6; MXD6; ROX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)