

## Product datasheet for **HP213596**

### COQ9 Human qPCR Primer Pair (NM\_020312)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57017
Forward Sequence:	ACTTGCGCTTCCTGGAAAACC
Reverse Sequence:	ACCCATGAGTCCTGCACCAGT
ACCN:	<u><a href="#">BC001478</a></u> , <u><a href="#">NM_020312</a></u> , <u><a href="#">NM_020312.1</a></u> , <u><a href="#">NM_020312.2</a></u> , <u><a href="#">NM_020312.3</a></u> , <u><a href="#">BC064946</a></u> , <u><a href="#">BC064946.1</a></u> , <u><a href="#">BC029341</a></u> , <u><a href="#">BC054340</a></u> , <u><a href="#">BM680448</a></u> , <u><a href="#">BM712215</a></u> , <u><a href="#">BM717991</a></u> , <u><a href="#">NM_020312.4</a></u>
UniProt ID:	<u><a href="#">O75208</a></u>
Synonyms:	C16orf49; COQ10D5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)