

Product datasheet for **HP213604**

VANGL2 Human qPCR Primer Pair (NM_020335)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57216
Forward Sequence:	AGCCTCAGTTCACGCTCAAGGT
Reverse Sequence:	GGTTGTAGACAGGGAAGTCATG
ACCN:	NM_020335 , NM_020335.1 , NM_020335.2 , BC103920 , BC166628 , BG059541 , BU162175 , BX420618 , NM_020335.3
UniProt ID:	Q9ULK5
Synonyms:	LPP1; LTAP; STB1; STBM; STBM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)