

Product datasheet for **HP213621**

PCNP Human qPCR Primer Pair (NM_020357)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57092
Forward Sequence:	CAGAGGAAATGCCTCCAGAAGC
Reverse Sequence:	ATTCGCTCCCACAGCTTCTGG
ACCN:	NM_020357 , NM_020357.1 , NM_020357.2 , BC022001 , BC022001.2 , BC013916 , BC107796 , BF244508 , BF541732 , BI547939 , BI596913 , BI913146 , BU608892 , BX647886 , NM_020357.3
UniProt ID:	Q8WW12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)