

## Product datasheet for **HP213623**

### phospholipid scramblase 2 (PLSCR2) Human qPCR Primer Pair (NM\_020359)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57047
Forward Sequence:	GTGGTCCATGTATCGTGTGCAG
Reverse Sequence:	AAGTTGTCAGCATCAGTAAATGCC
ACCN:	<u><a href="#">NM_020359</a></u> , <u><a href="#">NM_020359.1</a></u> , <u><a href="#">NM_020359.2</a></u> , <u><a href="#">BC069785</a></u> , <u><a href="#">BC055415</a></u> , <u><a href="#">BC120969</a></u> , <u><a href="#">BC141969</a></u> , <u><a href="#">BU852276</a></u>
UniProt ID:	<u><a href="#">Q9NRY7</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)