

## Product datasheet for **HP213647**

### HRAS like suppressor (HRASLS) Human qPCR Primer Pair (NM\_020386)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57110
Forward Sequence:	CGAAGTGTCCGTCCTGGCTAT
Reverse Sequence:	CAAGAGCTGCATTTTCACCAGGG
ACCN:	<a href="#">NM_020386</a> , <a href="#">NM_020386.1</a> , <a href="#">NM_020386.2</a> , <a href="#">NM_020386.3</a> , <a href="#">NM_020386.4</a> , <a href="#">BC048095</a> , <a href="#">BC048095.1</a> , <a href="#">BX098792</a> , <a href="#">NM_020386.5</a>
UniProt ID:	<a href="#">Q9HDD0</a>
Synonyms:	A-C1; H-REV107; HRASLS; HRASLS1; HRSL1; HSD28; PLA/AT1; PLAAT-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)