

Product datasheet for **HP213654**

BCL2L10 Human qPCR Primer Pair (NM_020396)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10017
Forward Sequence:	AGCACCGCGCCTGGCTGCAG
Reverse Sequence:	GACAGAAAAGCCTGGACCAGCT
ACCN:	<u>NM_020396</u> , <u>NM_020396.1</u> , <u>NM_020396.2</u> , <u>NM_020396.3</u> , <u>BC104442</u> , <u>BC093826</u> , <u>BC093828</u> , <u>BC104443</u> , <u>BC143227</u> , <u>BX119239</u>
UniProt ID:	<u>Q9HD36</u>
Synonyms:	BCL-B; bcl2-L-10; Boo; Diva
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)