

Product datasheet for **HP213656**

SPINLW1 (EPPIN) Human qPCR Primer Pair (NM_020398)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57119
Forward Sequence:	TATGGTGGCTGCCAGGAAACA
Reverse Sequence:	GCCACATTCTGGCAGTTCTTCC
ACCN:	<u>NM_020398</u> , <u>NM_020398.1</u> , <u>NM_020398.2</u> , <u>NM_020398.3</u> , <u>BC044829</u> , <u>BC044829.1</u> , <u>BC003675</u> , <u>BC053369</u> , <u>BF057305</u> , <u>NM_020398.4</u>
UniProt ID:	<u>O95925</u>
Synonyms:	CT71; CT72; dJ461P17.2; SPINLW1; WAP7; WFDC7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)