

Product datasheet for **HP213658**

NUP107 Human qPCR Primer Pair (NM_020401)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57122
Forward Sequence:	CACCATTGCCTGGAGTTGGCTA
Reverse Sequence:	GTAGTGCCAGTATCTAGGGCTG
ACCN:	NM_020401 , NM_020401.1 , NM_020401.2 , NM_020401.3 , BC043343 , BC043343.2 , BC017167 , BC017670 , BE326877 , NM_020401.4
UniProt ID:	P57740
Synonyms:	GAMOS7; NPHS11; NUP84; ODG6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)