

Product datasheet for **HP213661**

TEM1 (CD248) Human qPCR Primer Pair (NM_020404)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57124
Forward Sequence:	CTACGTTGGTGGCTTCGAGTGT
Reverse Sequence:	TCCAGCAACTCATCTCCGAGGT
ACCN:	NM_020404 , NM_020404.1 , NM_020404.2 , BC051340 , BC051340.1 , BC098158 , BC104484 , BC105633 , BQ005949 , NM_020404.3
UniProt ID:	Q9HCU0
Synonyms:	CD164L1; TEM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)