

Product datasheet for **HP213680**

SLC44A2 Human qPCR Primer Pair (NM_020428)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57153
Forward Sequence:	TTCTGGTGCTGGGCTACGGAAT
Reverse Sequence:	TCTGAAAGCCGAGGTCCACCAA
ACCN:	NM_020428 , NM_020428.1 , NM_020428.2 , NM_020428.3 , BC040556 , BC040556.1 , BC010617 , BC020511 , BC042692 , BI602883
UniProt ID:	Q8IWA5
Synonyms:	CTL2; PP1292
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)