

Product datasheet for **HP213689**

CaMKI gamma (CAMK1G) Human qPCR Primer Pair (NM_020439)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57172
Forward Sequence:	GGAGGACATCTATGAGAGCACC
Reverse Sequence:	TCACCAGACTGGCATCCTTCTC
ACCN:	<u>NM_020439</u> , <u>NM_020439.1</u> , <u>NM_020439.2</u> , <u>BC032787</u> , <u>BC032787.1</u> , <u>BG715920</u> , <u>NM_020439.3</u>
UniProt ID:	<u>Q96NX5</u>
Synonyms:	CLICK3; CLICKIII; dj272L16.1; VWS1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)